

Quala Dental Products

sharps



Made in USA

Sutures

Precision drilled end needle ensures smooth transition to suture material and maximizes secure attachment of suture to needle. Needles are sharpened to a smooth, precise cutting edge to penetrate cleanly and stay sharp for repeated passage through tissue. Double sterile pouch allows sterile transfer onto a sterile field. 12 per box.

FEATURES

Source for all dental sutures

Drilled-end needle

Unique needle design and manufacture

Double sterile pouch

See-through packaging

Industry standard label colors

BENEFITS

Can meet all dental suture requirements for single-source economy with the widest central suture range in the industry

Ensures smooth transition from needle to material; maximizes needle attachment strength

Needles stay sharp for repeated passage through tissue

Allows sterile transfer

Saves money by reducing selection errors and wasted materials

Simplifies selection, enhances efficiency



Quala
DENTAL PRODUCTS
www.quala.com



Quala Dental Products

sharps

ITEM#	NEEDLE CODE	NEEDLE DESCRIPTION	SUTURE SIZE	SUTURE LENGTH
Silk Black Braided				
P00066	C6	3/8 circle-19 mm reverse cutting	4/0	18" (45 cm)
P00072	C6	3/8 circle-19 mm reverse cutting	3/0	18" (45 cm)
P00062	C26	1/2 circle-15 mm reverse cutting	3/0	18" (45 cm)
P00064	C31	1/2 circle-24 mm reverse cutting	3/0	18" (45 cm)
Nylon Black Monofilament				
P00662	C6	3/8 circle-19 mm reverse cutting	4/0	18" (45 cm)
Plain Gut				
P00592	C6	3/8 circle-19 mm reverse cutting	3/0	18" (45 cm)
P00071	C6	3/8 circle-19 mm reverse cutting	4/0	18" (45 cm)
P00822	C6	3/8 circle-19 mm reverse cutting	3/0	27" (67.5 cm)
P00810	C7	3/8 circle-24 mm reverse cutting	3/0	27" (67.5 cm)
P00065	C31	1/2 circle-24 mm reverse cutting	3/0	18" (45 cm)
Chromic Gut				
P00687	C3	3/8 circle-19 mm reverse cutting	5/0	18" (45 cm)
P00558	C6	3/8 circle-19 mm reverse cutting	4/0	18" (45 cm)
P00525	C6	3/8 circle-19 mm reverse cutting	3/0	18" (45 cm)
P00635	C6	3/8 circle-19 mm reverse cutting	4/0	27" (67.5 cm)
P00636	C6	3/8 circle-19 mm	3/0	27" (67.5 cm)

Distributed by:



Quala
DENTAL PRODUCTS
www.quala.com

#QSS060R1

